Spirometer Software 12/4/2012

## designsafe Report

Application: Spirometer Software Analyst Name(s):

Description: Company:

Product Identifier: Facility Location:

Assessment Type: Detailed

Limits:

Sources:

Risk Scoring System: ANSI B11.0 (TR3) Two Factor

Guide sentence: When doing [task], the [user] could be injured by the [hazard] due to the [failure mode].

			Initial Assessment			Final Assessment		Status /
	User /	Hazard /	Severity		Risk Reduction Methods	Severity		Responsible
Item Id	Task	Failure Mode	Probability	Risk Level	/Comments	Probability	Risk Level	/Reference
1-1-1	All Users	mechanical: cutting	Minor	Low	making the mouthpiece with	Minor	Negligible	
	Common Tasks	sharp mouthpiece	Likely		smooth edges	Unlikely		
1-1-2	All Users Common Tasks	mechanical : impaired movement	Minor Likely	Low	warning label(s)	Minor Likely	Low	
		tangling in cord	,					
1-1-3	All Users	mechanical : swallowing	Moderate	Low	strong connection between	Moderate	Low	
	Common Tasks	foreign object swallowing microphone, swallowing foreign object in	Unlikely		components, warning label(s), clean between users	Unlikely		
		microphone cavity						
1-1-4	All Users Common Tasks	mechanical : stangulation tangling in cord	Catastrophic Unlikely	Medium	warning label(s)	Catastrophic Unlikely	Medium	
1-1-5	All Users Common Tasks	mechanical : chocking hazard inhaling microphone, inhaling foreign object in microphone cavity	Minor Unlikely	Negligible	strong connection between components, warning label(s), clean between users	Minor Unlikely	Negligible	
1-1-6	All Users Common Tasks	electrical / electronic : electrocution exposed wires	Catastrophic Unlikely	Medium	fixed enclosures / barriers	Catastrophic Remote	Low	
1-1-7	All Users	slips / trips / falls : tripping	Minor	Medium	warning label(s)	Minor	Low	
	Common Tasks	loose cord	Very Likely			Likely		
1-1-8	All Users Common Tasks	slips / trips / falls : falling spirometer dropping spirometer	Minor Very Likely	Medium	substitute less hazardous material / methods	Minor Very Likely	Medium	

Spirometer Software 12/4/2012

Item Id	User / Task	Hazard / Failure Mode	Initial Assessn Severity Probability	nent Risk Level	Risk Reduction Methods /Comments	Final Assessm Severity Probability	ent Risk Level	Status / Responsible /Reference
1-1-9	All Users Common Tasks	slips / trips / falls : falling computer dropping computer, pulling computer off surface through cord	Minor Likely	Low	make connection cord that has adjustable length	Minor Likely	Low	
1-1-10	All Users Common Tasks	ergonomics / human factors : excessive force / exertion too many tests	Minor Likely	Low	limit the number of tests in a certain time period to be performed by a single user	Minor Unlikely	Negligible	
1-1-11	All Users Common Tasks	ergonomics / human factors : repetition too many tests	Minor Likely	Low	limit the number of tests in a certain time period to be performed by a single user	Minor Unlikely	Negligible	
1-1-12	All Users Common Tasks	environmental / industrial hygiene : irritants foreign objects in spirometer cavity	Minor Likely	Low	clean between users	Minor Unlikely	Negligible	
1-1-13	All Users Common Tasks	environmental / industrial hygiene : cross-contamination sanitation failure	Moderate Likely	Medium	instruct to wash hardware between users, place disposable filter on microphone	Moderate Remote	Negligible	
1-1-14	All Users Common Tasks	Other : emotional stress false positives, false negatives	Moderate Likely	Medium	add in accuracy and precision measures, ensure that only doctors are using to diagnose, instruction manuals	Moderate Unlikely	Low	
1-1-15	All Users Common Tasks	Other : security of patient information security breach	Serious Unlikely	Medium	compliane with HIPAA regulations	Serious Unlikely	Medium	
1-1-16	All Users Common Tasks	Other: software error overburdened computer, faulty sound file	Minor Unlikely	Negligible	visible alarm or signal	Serious Unlikely	Medium	
1-1-17	All Users Common Tasks	Other: overburdened user computer too large of a downloadable component	Minor Unlikely	Negligible	visible alarm or signal	Minor Remote	Negligible	

Spirometer Software 12/4/2012

Item Id	User / Task	Hazard / Failure Mode	Initial Assessmer Severity Probability	nt Risk Level	Risk Reduction Methods /Comments	Final Assessmen Severity Probability	t Risk Level	Status / Responsible /Reference
1-1-18	All Users Common Tasks	Other: incorrect test execution misunderstanding of device operation, lack of proper guidance	Moderate Likely	Medium	instruction manuals	Moderate Unlikely	Low	
1-1-19	All Users Common Tasks	Other : wrong diagnosis misunderstanding of results	Moderate Likely	Medium	instruction manuals	Moderate Unlikely	Low	